

Setting up CODE CR1500 Scanner for Rapid Tags as a Virtual COM Port

The instructions below will help you download and install the CODE Virtual COM Port Drivers and program your CODE Scanner. More information may be found on the Rapid Lab website: <http://rapidlab.com/scanner-and-software-support/>

1. Install the correct USB Virtual COM Port driver for your system:

- a. Visit the following link to see all products.
<https://ww2.codecorp.com/downloads.php>
- b. Click on the correct product to view available firmware and software downloads.
- c. Scroll to the bottom of the page and check the agreement box for the correct firmware or driver and click on the download button.
- d. **Once the driver has installed, restart your computer.**

2. Reset the Scanner:

- a. **Unplug** the CODE CR1500 from the PC.
- b. **While holding the trigger** pressed, plug the reader into the USB port on the PC and hold for 5-10 seconds. You will hear 5 rapid beeps and then you can let go of the trigger. The unit is now in Boot Mode.
- c. Next, **hold the trigger down again** for approximately 5-10 seconds until it reboots. The reader will vibrate first and then beep several times and is now reset.

3. Reprogram the scanner:

- a. **Scan** the codes below to reprogram your scanner.

Required Codes:

Custom Settings



CC004482_03

Enable VCOM and Disable Data Formatting



CC004854



Additional optional codes on the following page below.

Optional Codes:

1. Continuous Scan On



M20224_01

2. 1 Second Duplicate Scan Delay



M20230_01

3. 3 Second Duplicate Scan Delay



M20232_01

4. Disable Duplicate Scan Delay



M20229_01

1. This code allows for continuous scanning on your CR1500.

2. This code creates a 1 second delay between scans.

3. This code creates a 3 second delay between scans.

4. This code disables duplicate scanning.