

## Setting up CODE CR1500 Scanner for Rapid Tags as a Virtual COM Port

The instructions below will help you download and install the CODE Virtual COM Port Drivers and program your CODE Scanner. More information may be found on the Rapid Lab website: <http://rapidlab.com/scanner-and-software-support/>

### 1. Install the correct USB Virtual COM Port driver for your system:



- Visit the following link and click on the 'Software' navigation tab.  
<https://codecorp.com/products/code-reader-1500>
- Scroll down and click on 'USB Virtual COM Driver (32-bit and 64-bit Operating Systems)' then click on Download.
- Download and install the driver.
- Once the driver has installed, restart your computer.**

### 2. Reset the Scanner:

- Unplug** the CODE CR1500 from the PC.
- While holding the trigger** pressed, plug the reader into the USB port on the PC and hold for 5-10 seconds. You will hear 5 rapid beeps and then you can let go of the trigger. The unit is now in Boot Mode.
- Next, **hold the trigger down again** for approximately 5-10 seconds until it reboots. The reader will vibrate first and then beep several times and is now reset.

### 3. Reprogram the scanner:

- Scan** the codes below to reprogram your scanner.

Required Codes:	
<p>Custom Settings</p>  <p>CC004482</p>	<p>Enable VCOM and Disable Data Formatting</p>  <p>CC004854</p>